



**OraSure Technologies, Inc.**

diagnostic solutions for the new millennium

## **OraQuick® Rapid HIV-1 Antibody Test Executive Summary**

### **Product Description:**

The OraQuick® Rapid HIV-1 Antibody Test is a single-use, qualitative immunoassay to detect antibodies to Human Immunodeficiency Virus Type 1 (HIV-1) in fingerstick whole blood specimens. The OraQuick® Rapid HIV-1 Antibody Test is intended for use as a point-of-care test to aid in the diagnosis of infection with HIV-1. This test is suitable for use in multi-test algorithms designed for statistical validation of rapid HIV test results. When multiple rapid HIV tests are available, this test should be used in appropriate multi-test algorithms.

### **Price:**

\$12.00 U.S. This introductory price is valid through February 28, 2003, after which the price per kit goes to our list price of \$20.00. Shipping is paid by the customer and disclosed once a ship to address and zip code are provided.

### **Tax Identification Number:**

# 36-4370966

### **National Stock Number:**

100 test kit is: 6550-01-502-6235

25 test kit is: 6550-01-502-6241

### **Material Provided:**

The OraQuick® Rapid HIV-1 Antibody Test, 25 count, test kit - includes 25 devices, 30 loops, 6 test stands, 1 product instruction, 25 subject information pamphlets, and 1 customer letter. The dimensions and weight are (25 Count Box) 12 x 9 x 6 and 2 lbs, respectively.

### **Materials Required and Available As An Accessory To The Kit:**

OraQuick® Rapid HIV-1 Antibody Test Kit Controls. Test Kit Controls contain Positive Control (1 vial, black cap, 0.2 mL) and Negative Control (1 vial, white cap, 0.2 mL), and a Package Insert

### **Material Required, but not provided:**

Timer or watch  
Antiseptic wipe  
Sterile lancet  
Sterile gauze pads  
Disposable gloves  
Biohazard disposal container

### **Storage:**

Store unused OraQuick® Rapid HIV-1 Antibody Tests unopened at 2°–27°C (35–80°F). Do not open the Divided Pouch until you are ready to perform a test. If stored refrigerated, ensure that the Divided Pouch is brought to ambient temperature (15°–27°C) before opening.

### **Shelf Life:**

7 months

David C. Oxley

<mailto:dcoxley@orasure.com>

150 Webster Street Bethlehem, Pennsylvania 18015-1389

Phone 484-695-8570 Fax 610-882-2275

[www.orasure.com](http://www.orasure.com)